IVD

IRMA IGF-I

REF

A15729

**IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION
OF INSULIN-LIKE GROWTH FACTOR I IN HUMAN SERUM AND PLASMA**

1. PRINCIPLE OF THE ASSAY

The immunoradiometric assay of Insulin-like Growth Factor I (IGF-I) is a "sandwich" type assay. Mouse monoclonal antibodies directed against two different epitopes of IGF-I and hence not competing are used. In order to release IGF-I from its binding proteins, a prior dissociation step is necessary. Samples and calibrators are incubated in tubes coated with the first monoclonal antibody in the presence of the second monoclonal antibody labeled with iodine 125. After incubation, the content of tubes is removed and bound radioactivity is measured. Unknown values are determined by interpolation from a standard curve. The bound radioactivity is directly proportional to the IGF-I concentration in the sample.

2. DESCRIPTION OF REAGENTS

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Storage conditions for reagents after reconstitution or dilution are indicated in paragraph Assay Procedure.

2.1 Anti-IGF-I antibody-coated tubes : 2 x 50 tubes (ready-to-use)

2.2 ¹²⁵I-labeled anti-IGF-I antibody : one 33 mL vial (ready-to-use)

The vial contains 370 kBq, at the date of manufacture, of ¹²⁵I-labeled immunoglobulins in liquid form, containing bovine serum albumin, sodium azide (<0.1%, see § Precautions), and a dye.

2.3 Calibrators : six vials (lyophilized)

The calibrator vials contain from 0 to 1,200 ng/mL of IGF-I in a buffer with bovine serum albumin and a preservative. The exact concentration is indicated on each vial label. The calibrators were calibrated against the international reference standard, WHO 87/518.

2.4 Control sample : one vial (lyophilized)

The vial contains IGF-I lyophilized in bovine serum albumin. The expected values are in the concentration range indicated on the vial label.

2.5 Dissociation buffer : two 25 mL vials (ready-to-use)

The vial contains bovine serum albumin.

2.6 Wash solution (20x) : one 50 mL vial

Concentrated solution has to be diluted before use.

Note : Temperatures and expiry dates printed on vial labels only apply to the long-term storage of components by the manufacturer. Do not take into account.

3. MATERIAL REQUIRED BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- Precision micropipet (50 µL).
- Semi-automatic pipets (300 µL, 1 mL).
- Vortex type mixer.
- Horizontal or orbital shaker (min. 300 rpm).
- Aspiration system.
- Gamma counter set for 125 iodine.

4. PRECAUTIONS

4.1 General remarks

- Bring all reagents to room temperature before pipeting.
- Do not mix the reagents from kits of different lots.
- A standard curve must be included with each assay.
- The correct setting of the shaker is very important for the reproducibility of the assay.
- It is recommended to perform the assay in duplicate.
- Strictly observe the order in which reagents are to be added.
- Each tube must be used only once.

4.2 Protection against ionizing radiation

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use.

Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

4.3 Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

4.4 Human serum

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS. Waste should be discarded according to the country rules.

5. SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in dry tubes or in tubes with heparin or EDTA.
- Separate serum or plasma from cells by centrifugation.
- Serum and plasma samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen (<-18°C) after aliquoting so as to avoid repeated freezing and thawing.
- If samples have concentrations greater than the highest calibrator, they must be diluted with the dissociation buffer.

6. ASSAY PROCEDURE

6.1 Preparation of reagents

Let all the reagents come to room temperature.

6.1.1 Reconstitution of calibrators and control sample

The contents of the vials must be brought to room temperature before reconstitution with the volume of distilled water indicated on the vial label. Wait for 30 min following reconstitution and mix gently to avoid foaming before dispensing. Store the reconstituted solutions at 2-8°C for one week or aliquoted at <-18°C for a longer time, until the expiry date of the kit.

6.1.2 Preparation of the wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution may be stored at 2-8°C until the expiry date of the kit.

6.2 Treatment of samples and control

Do not treat calibrators.

- To plastic tubes, add successively 50 µL of sample and 1 mL of dissociation buffer (these conditions allow to treat 50 samples; divide volumes by two to treat 100 samples, i.e. 25 µL of serum and 500 µL of buffer).
- Mix with a vortex-type shaker.
- Treated samples may be kept for 48h at 2-8°C; for longer storage keep at <-18°C.

6.3 Assay procedure (see table next page).

7. RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of IGF-I concentrations in samples measured at the same time as the calibrators.

7.1 Standard curve

The results in the package insert were calculated using a semi-logarithmic curve fit ("spline" mode) with B/T (%) or B/Bmax (%) on vertical axis and the IGF-I concentration of the calibrators on the horizontal axis (ng/mL). Other data reduction methods may give slightly different results.

Total activity : 127,559 cpm				
Calibrators	IGF-I ng/mL	cpm (n=3)	B/T (%)	B/Bmax (%)
0	0	257	0.20	0.43
1	30	1,138	0.89	1.93
2	100	4,496	3.52	7.65
3	300	18,620	14.59	31.72
4	600	39,285	30.79	66.93
5	1,200	58,695	46.00	100

(Example of standard curve, do not use for calculation)

7.2 Samples

Locate the ratio B/T (%) or B/Bmax (%) on the vertical axis of the standard curve, read-off the IGF-I concentration of the sample on the horizontal axis in ng/mL.

8. QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended that their results be analyzed using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address:

immuno-techsup@beckmancoulter.com

9. EXPECTED VALUES

It is suggested that each laboratory establish its own normal values. The following values obtained with healthy subjects are indicative only.

Adults

Age (years)	n	IGF-I (ng/mL)				
		Min.	5 th percentile	Median	95 th percentile	Max.
20-30	51	219	232	288	385	644
30-40	44	140	177	245	382	405
40-50	43	64	124	199	290	336
50-60	18	71	71	147	263	284
60-70	20	94	94	141	269	269
70-80	20	72	76	117	160	167

Children

Stage of puberty (years)	n	IGF-I (ng/mL)			
		Mean	Min.	Max.	
P1	0-4	5	114	49	171
	>4	27	250	76	499
P2	6	303	247	396	
P3	7	414	249	642	
P4-P5	7	400	271	550	

Children constitutionally small

These children have a height that is lower by two standard deviations or more than average with GH concentration greater than 20 mIU/L after stimulation, and a regular rate of growth.

Stage of puberty (years)	n	IGF-I (ng/mL)			
		Mean	Min.	Max.	
P1	0-4	13	114	98	180
	5-7	25	115	98	156
	8-9	21	129	76	186
	10-11	24	151	76	234
	>12	27	198	131	278
P2	20	258	163	502	
P3	14	351	185	617	
P4	9	423	272	540	

10. PERFORMANCE CHARACTERISTICS

(for more details, see the data sheet "APPENDIX")

10.1 Analytical sensitivity: 2 ng/mL

10.2 Specificity

The antibodies used in the immunoassay are highly specific for IGF-I. Extremely low cross reactivities were obtained against several molecules (insulin, proinsulin, IGF II, GH).

10.3 Precision

10.3.1 Intra-assay

Samples were assayed in 25 times in the same series. The coefficients of variation were found below or equal to 6.3%.

10.3.2 Inter-assay

Samples were assayed in duplicate in 25 different series. Coefficients of variation were found below or equal to 6.8%.

10.4 Accuracy

10.4.1 Dilution test

High-concentration samples were serially diluted with the zero calibrator. The recovery percentages obtained were between 88% and 111%.

10.4.2 Recovery test

Low-concentration samples were spiked with known quantities of IGF-I. The recovery percentages obtained were between 91% and 103%.

10.5 Measurement range (from analytical sensitivity to highest calibrator):

2 – 1,200 ng/mL.

11. LIMITATIONS OF THE PROCEDURE

The non-respect of the instructions in this package insert may affect results significantly. Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

6.3 ASSAY PROCEDURE

Step 1 Additions	Step 2 Incubation	Step 3 Counting
<p>To antibody coated tubes, add successively in this order:</p> <ul style="list-style-type: none"> - 300 µL of tracer*, - 50 µL of calibrator, control or sample, <p>Mix.</p>	<p>Incubate 60 minutes at 18-25°C with shaking (350 rpm).</p>	<p>Aspirate carefully the content of tubes (except the 2 tubes «total cpm»).</p> <p>Wash twice with 2 mL of wash solution.</p> <p>Count bound cpm (B) and total cpm (T) for 1 min.</p>

* Add 300 µL of tracer to 2 additional tubes to obtain total cpm.

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REF

A15729

IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF
INSULIN-LIKE GROWTH FACTOR I IN HUMAN SERUM AND PLASMA

Performance characteristics

A. Specificity

The specificity of the assay was determined by measuring the apparent IGF-I value given by high concentrations of related compounds in the absence (cross-reactivity) or presence (interference) of IGF-I.

A.1 Cross-reactivity

Related molecules	Concentration	IGF-I concentration measured	
	nM	ng/mL	nM
GH	40	0	0
Insulin	20	0	0
Proinsulin	40	0	0
IGF-II	650	0	0

A.2 Interferences

Spiking	IGF-I concentration expected		IGF-I concentration measured	
	ng/mL	nM	ng/mL	nM
GH (40 nM) + IGF-I	270	35.3	274	35.6
Insulin (25 nM) + IGF-I	270	35.3	273	35.6
Proinsulin (20 nM) + IGF-I	270	35.3	266	34.8
IGF-II (650 nM) + IGF-I	294	38.5	278	36.1

B. Accuracy

B.1 Dilution test

Three samples containing IGF-I were treated with the dissociation buffer following the procedure indicated in the insert. The treated samples were then diluted using the dissociation buffer.

Sample	Dilution factor	IGF-I (ng/mL)		Recovery (%) measured/expected
		measured	expected	
1	Undiluted	917	-	-
	1/2	415	459	90.4
	1/4	203	229	88.6
	1/8	113	115	98.3
	1/16	57	57	100
2	Undiluted	918	-	-
	1/2	440	459	95.8
	1/4	223	230	96.9
	1/8	115	115	100
	1/16	60	57	105.3
3	Undiluted	779	-	-
	1/2	386	390	99.0
	1/4	192	195	98.5
	1/8	99	97	102.1
	1/16	54	49	110.2

B.2 Recovery test

IGF-I was added to serum samples and assayed according to the procedure of the kit.

Sample (ng/mL)	IGF-I (ng/mL)		Recovery measured/expected (%)
	measured	expected	
1 (61)	530	561	94.5
	313	311	100.6
	192	186	103.2
	139	136	102.2
2 (109)	557	609	91.5
	347	359	96.7
	224	234	95.7
	180	184	97.8
3 (79)	525	579	90.7
	314	329	95.4
	195	202	96.5
	157	154	101.9

C. Precision

Samples	Intra-assay (n=25)				Inter-assays (n=25)			
	A	B	C	D	E	F	G	H
Mean value (ng/mL)	40	81	229	575	25	82	129	244
CV (%)	3.2	6.3	2.4	3.8	6.8	5.7	5.3	6.8

D. Drift

The choice of optimal conditions for the immunoassay has minimized drift. A 30-minutes lag between the addition of the tracer and addition of the samples does not affect the assay.

E. Hook effect

No hook effect was observed up to a concentration of IGF-I of approximately 7,000 ng/mL.

F. Effect of anti-coagulants

The effect of differences in collection of blood were tested on samples of normal individuals. The mode of collecting sample does not affect assay.

Samples	IGF-I ng/mL							Ratio Dry tubes / others
	1	2	3	4	5	6	7	
Dry	191	239	213	150	269	203	146	100%
with gel	185	245	216	167	272	207	144	101%
with heparin	185	239	183	200	287	204	138	101%
with EDTA	184	229	163	179	256	204	135	95%
with Citrate	159	223	149	160	220	183	125	87%

$T_{1/2} (^{125}\text{I}) = 1443 \text{ h} = 60.14 \text{ d}$

¹²⁵ I	E (MeV)	%
γ	0.035	
X	0.027	114
	0.032	25

Use by / Utiliser jusqu'à / Verwendbar bis / Utilizzare entro / Use antes de / Ημερομηνία λήξης / Použitelné do / Použitelné do / Užyc do / Felhasználható / Naudoti iki / Использовать до

IVD

In Vitro diagnostic Device / Diagnostic in Vitro / In Vitro Diagnostikum / Diagnostico in vitro / Diagnostico in Vitro / In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν / Pro diagnostiku in vitro / Pre diagnostiku in vitro / Do diagnostyki in vitro / In vitro diagnosztikum / In vitro diagnostikai / Для ин vitro диагностики

REF

Catalogue Number / Référence catalogue / Bestellnummer / Numero di catalogo / Numero de catalogo / Αριθμός καταλόγου / Katalogové číslo / Katalógové číslo / Numer katalogowy / Katalógusszám / Katalogo Nr. / Καταложный №

LOT

Batch code / Numéro de lot / Chargenbezeichnung / Coddice del lotto / Código de lote / Αριθμός Παρτίδας / Číslo šarže / Číslo šarže / Numer serii / Lot szám / Serijos Nr. / № серии

Caution, see instructions for use / Attention, consulter la notice d'utilisation / Achtung, Gebrauchsanweisung beachten / Attenzione vedere le istruzioni per l'uso / Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα / Pozor, čtěte pozorně návod / Pozor, čítajte návod rozorne / Uwaga, patrz instrukcja użycia / Figyelem, olvassa el a használati utasítást / Dėmesio, atidžiai perskaitykite instrukciją / Внимание, читайте тщательно инструкцию

Consult Instructions for Use / Consulter la notice d'utilisation / Gebrauchsanweisung beachten / Leggere le istruzioni per l'uso / Consulte las instrucciones de uso / Συμβουλευτείτε τις οδηγίες χρήσης / Viz návod / Vid návod / Patrz instrukcja użycia / Olvassa el a használati utasítást / Žiūrėti instrukciją / Смотри инструкцию

Temperature limitation / Limites de température / Temperaturgrenzen / Limiti di temperatura / Limites de temperatura / Περιορισμοί θερμοκρασίας / Rozmezi skladovacích teplot / Rozsahy skladovacích teplôt / Ograniczenie temperatury / Hőmérséklethatárok / Saugojimo temperatūrų intervalas / Диапазон температур хранения

Radioactive / Radioactif / Radioattivo / Radioaktiv / Radioactivo / Ραδιενεργός / Radioaktivní / Rádioaktívne / Radioaktywność / Radioaktív / Radioaktyvus / Радиоактивный

For XX tests / Pour XX dosages / Für XX Bestimmungen / Per XX dosaggi / Para XX ensayos / Περιεχόμενο επαρκές για XX εξετάσεις / Lze použít pro XX testů / Určené na XX testov / Zawartość na XX testów / Tartalma XX teszt elvégzésére elegendő / Pakanka XX testui atlikti / Предназначен для XX тестов

Manufacturer / Fabricant / Hersteller / Fabbicante / Fabricante / Κατασκευαστής / Výrobce / Výrobca / Producent / Gyártó / Gamintojas / Производитель