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

Immune-mediated peripheral neuropathies	Number of samples	Anti-ganglioside autoantibody profiles detected by in-house immunodot	Number of positive samples using minimal positivity criteria*				
			Ih T	T1	T2	T3	T4
AMAN and AMSAN	6	IgG antibodies against GM1 > GD1b	6	6	3	1	6
Miller Fisher syndrome	4	IgG antibodies against GQ1b	4	4	2	3	3
GBS with ophthalmoplegia	1	IgG antibodies against GD1a, GQ1b	1	1	0	0	0
GBS with ophthalmoplegia	1	IgG antibodies against GD1a, GD1b, GT1b, GQ1b	1	1	0	0	1
GBS post CMV infection	2	IgM antibodies against GM2	2	2	2	2	2
Axonal motor GBS post CJ infection	1	IgG and IgM antibodies against GM1 and GD1b	1	1	0	0	1
Total of acute PN	15		15	15	7	6	13
CANOMAD syndrome	12	M-IgM antibodies against GD1b, GD3, GT1b, GQ1b	12	12	6	5	7
Chronic motor neuropathy	3	M-IgM antibodies against GM1 and GD1b	3	3	2	1	3
Chronic motor neuropathy with lymphoma	1	M-IgM antibodies against GD1a and GT1b	1	1	1	1	1
Chronic sensory neuropathy	1	M-IgM antibodies against sulfatides > GD1b > GM1	1	1	1	1	1
MMN	1	IgM antibodies against GM1 > GD1b	1	1	0	0	1
Total of chronic PN	18		18	18	10	8	13
Total of PN	33		33	33	17	14	26

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1. Conrad K, Schöbeler W, Hiepe F, Fritzler MJ Autoantibodies in organ specific autoimmune diseases - a diagnostic reference. Autoantigens, Autoantibodies, Autoimmunity 2011, Volume 8
2. Yuki N & Hartung HP. Guillain-Barré syndrome. N Engl J Med. 2012;366(24):2294-304.
3. Ang CW, Jacobs BC, Laman JD. The Guillain-Barré syndrome: a true case of molecular mimicry. Trends Immunol. 2004;25(2):61-6.
4. Willison HJ & Yuki N. Peripheral neuropathies and anti-glycolipid antibodies. Brain. 2002;125(Pt 12):2591-625.
5. Caudie C, Quittard Pinon A, Bouhour F, Vial C, Garnier L, Fabien N. Comparison of commercial tests for detecting multiple anti-ganglioside autoantibodies in patients with well-characterized immune-mediated peripheral neuropathies. Clin Lab. 2013;59(11-12):1277-87.

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